

The influence of storage time on ponazuril concentrations in feline plasma

Sherry Cox ^{Corresp., 1}, Lainey Harvill ¹, Sarah Singleton ¹, Joan Bergman ¹, Becky DeBolt ²

¹ Department of Biomedical and Diagnostic Sciences, College of Veterinary Medicine, University of Tennessee, Knoxville, TN, United States

² Department of Small Animal Clinical Sciences, College of Veterinary Medicine, University of Tennessee, Knoxville, TN, United States

Corresponding Author: Sherry Cox

Email address: scox6@utk.edu

Background. The pharmacokinetics of ponazuril have been determined in several species however, there is very little information on the stability of the drug after storage for long periods of time. This study was undertaken to determine the stability of ponazuril in plasma samples stored at -80°C, which is the temperature most commonly used in the author's laboratory.

Method. Spiked plasma samples (0.3, 7.5, and 15 µg/mL) were stored at -80°C for three months. Analysis occurred on the first day and then once a week for the following twelve weeks. The drug was extracted using a chloroform extraction and separated by high performance liquid chromatography using ultraviolet detection.

Results. There was no loss of drug for any concentration for the first four weeks of storage. There was an average loss of less than 5% from day 35 through day 70 and an average loss of 6% on day 77 and 84. The data suggest that ponazuril is stable for 4 weeks when stored at -80°C and undergoes minimal loss in the remaining 8 weeks.

The Influence of Storage Time on Ponazuril Concentrations in Feline Plasma

Sherry Cox^{*1}, Lainey Harvill¹, Sarah Singleton¹, Joan Bergman¹, Becky DeBolt²

¹Department of Biomedical and Diagnostic Sciences and

²Department of Small Animal Clinical Sciences University of Tennessee Knoxville, College of Veterinary Medicine, Knoxville, TN, USA

Corresponding Author:

Sherry Cox, 2407 River Drive, Knoxville, TN, 37996 USA

Email address: scox6@utk.edu

Abstract

Background. The pharmacokinetics of ponazuril have been determined in several species however, there is very little information on the stability of the drug after storage for long periods of time. This study was undertaken to determine the stability of ponazuril in plasma samples stored at -80°C, which is the temperature most commonly used in the author's laboratory.

Method. Spiked plasma samples (0.3, 7.5, and 15 µg/mL) were stored at -80°C for three months. Analysis occurred on the first day and then once a week for the following twelve weeks. The drug was extracted using a chloroform extraction and separated by high performance liquid chromatography using ultraviolet detection.

Results. There was no loss of drug for any concentration for the first four weeks of storage. There was an average loss of less than 5% from day 35 through day 70 and an

average loss of 6% on day 77 and 84. The data suggest that ponazuril is stable for 4 weeks when stored at -80°C and undergoes minimal loss in the remaining 8 weeks.

Introduction

Ponazuril (also known as toltrazuril sulfone) is from the triazine class of antiprotozoal drugs and is a metabolite of the drug toltrazuril (Love, et al., 2015). Triazine compounds have been used in veterinary medicine to treat the harmful effects of protozoa such as *Toxoplasma gondii*, *Coccidia* spp., and *Isopora*, (Stock, et al., 2018).

The use of drugs in different species requires knowledge of their pharmacokinetics. An understanding of the pharmacokinetics is necessary for determining the most appropriate dose, route of administration and dosing interval. The determination of accurate drug concentrations is extremely important in the estimation of pharmacokinetic parameters and may impact design of the drug regimen. To determine these properties, it is essential to calculate the concentration of the compound of interest after administration. Because of the importance placed on the validation of bioanalytical methods that are used in the determination of data for pharmacokinetic studies the FDA has established guidelines (FDA, Guidance for industry: Bioanalytical Method Validation). The information in the guidelines applies to bioanalytical procedures such as chromatographic methods. A necessity for the generation of the data is a validated analytical procedure. However, if proper conditions are not maintained for storage the results may not accurately indicate the sample concentration at the time it was determined. For most studies, samples are stored until there is an adequate quantity for analysis, which requires them to be stored, usually

frozen for a period of time. To guarantee accurate results it is essential to determine drug degradation from samples using the storage conditions anticipated for the study.

Pharmacokinetic studies of ponazuril have been determined in several species (Laczay, et al., 1995; Epe, et al., 2005; Ghanem, et al., 2008; Love, et al., 2015; Prado, et al., 2011; Dirikolu, et al., 2011; Zou, et al., 2014; Stock, et al., 2018; Furr and Kennedy, 2020). To the author's knowledge there are no other stability articles involving drugs used in felines or triazine drugs. Thus, the intention of this research was to determine if concentrations of ponazuril in feline plasma decreased after storage for 12 weeks at -80°C using a validated method (Cox, et al., 2021).

Materials & Methods

Equipment

The equipment used for ponazuril assessment consisted of a 2695 separation module and a 2487 ultraviolet detector (Waters Milford, MA, USA). Separation was achieved on a Waters Symmetry Shield RP₁₈ (4.6 mm x 150 mm, 5 µm) column. The mobile phase consisted of 0.1% formic acid in water and acetonitrile (50:50 v/v). The ultraviolet detector was set at 254 nm and the flow rate was 1.1 mL/min.

Reagent and solutions

Ponazuril and diclazuril were purchased from Cayman Chemical (Ann Arbor, MI, USA). Both compounds had a purity greater than 99%. All other solvents and chemicals were acquired from Fisher Scientific (Pittsburg, PA, USA). Stock solutions of ponazuril (100 µg/mL) and diclazuril (100 µg/mL), which is the internal standard, were prepared by dissolving appropriate amounts of each compound in methanol. Cat plasma was

fortified with ponazuril to construct daily standard curves (0.1 – 25 µg/mL). These were accepted if the coefficient of determination (r^2) was greater than 0.99.

Sample Collection

Six adult cats from UTCVM faculty and staff were determined to be healthy based on results of physical examination and history were used. Healthy animals are routinely used in pharmacokinetic studies in order to prevent distortion of results due to abnormal metabolism, disease processes or interference from drugs the animal is taking. Venous blood was collected from the jugular vein of each cat and placed in lithium heparin tubes. The volume of blood collected from each cat did not exceed 1% of body weight in kilograms. Enough collected blood was centrifuged to provide plasma for the study. Plasma was pooled into one container and then divided into 3 separate tubes, and then spiked with one of the ponazuril calibration standards (0.3, 7.5, and 15 µg/mL) that were used in the study. Aliquots from each of the 3 individual standards were then placed in separate vials labeled Day 1, 7, 14, 21, 28, 35, 42, 49, 56, 63, 70, 77, and 84 for analysis on those dates. These concentrations were chosen as low medium and high standards within the linear curve used. Samples were then placed in the -80°C. Day 1 samples were analyzed on that day. Because of animal size and availability enough plasma was obtained to analyze only one sample for each ponazuril concentration evaluated.

Extraction method

Ponazuril was extracted from plasma samples that were thawed one time using a liquid-liquid extraction method (Cox et al., 2021). Frozen plasma samples were thawed at room temperature (22°C) for ~30 minutes, vortex-mixed, and 100 µl of plasma was

transferred to a screw top tube followed by 10 µl of diclazuril (100 µg/mL, internal standard) and 2 mL chloroform. The tubes were placed on a rocker for 15 minutes and then underwent centrifugation for 20 min at 1000 x g. The organic layer was transferred to a glass tube and evaporated to dryness with nitrogen gas. Samples were reconstituted in 250 µL of mobile phase and 100 µL was analyzed.

Statistical analysis

Comparative changes in the actual values of different measured ponazuril concentrations in plasma samples were used to summarize the effect of storage length on stability of the drug.

Results

Samples were analyzed on days 1, 7, 14, 21, 28, 35, 42, 49, 56, 63, 70, 77, and 84 (Figure 1A-D and Table 1). The actual values of the different measured ponazuril concentrations and the percent variation of the concentration compared to the starting concentration are reported in Table 1. Ponazuril concentrations were determined in all samples. There was no loss of drug for any concentration (low, medium or high) for the first four weeks (days 1 through 28) of storage at -80°C. Please refer to the Figure and Table for the remainder of the results.

The ponazuril method of analysis that was used for the stability study was previously validated (Cox et al., 2021). The method produced a linear curve over the concentration range used (0.1 – 25 µg/ml) with a lower limit of quantification of 0.1 µg/mL. The intraday and interday variability ranged from 3.7 to 10% and 2.1 to 5.5% while the average ponazuril recovery was 99% and the average recovery of diclazuril was 94%. Calibration curves were analyzed each day of plasma analysis.

Discussion

In order to determine the optimal storage time for ponazuril samples at -80°C, plasma was collected from healthy felines and pooled in order to reduce the effects of inter-individual variability. A validated method of analysis was applied to the samples to determine the effect of storage duration. Plasma samples spiked with various amounts of ponazuril that fell within a validated linear concentration range and stored in a -80°C freezer for 84 days were used. One limitation to the study was a lack of duplicate samples for the analysis. The analysis of only one sample did not allow for any statistical analysis.

Ponazuril was detected in all concentrations for the entire 84 days. For the first 28 days, there was no change in any of the concentrations. Indicating that the drug is stable after storage for that length of time. After 28 days, there was a slight loss of drug for all three concentrations. The average loss was 4% or less for days 35 through 63, which would still not have an impact on higher concentrations and have a minimal impact on lower concentrations. The samples should still be viable for this period of time. On day 70, there was an average loss of 5% and on days 77 and 84 there was an average loss of 6%. The loss on day 84 was greater for the 0.3 and 7.5 than the 15 µg/mL standard. Overall, the average loss from day 7 though day 84 was 2%, 4% and 2% for 0.3, 7.5 and 15 µg/mL, respectively.

To the author's knowledge there are no other stability articles involving drugs used in felines or triazine drugs. There is a paucity of information in the literature about the storage stability of drugs in biological samples in general. Determination of the concentrations in biological samples is fundamental to describing pharmacokinetic

parameters which in turn leads to appropriate dosing regimens. Therefore, it is imperative to determine the effect that storage can have on drug stability.

Conclusions

The intention of this study was to determine the effect of storage duration at -80°C on the stability of ponazuril concentrations in plasma. In summary, plasma samples, containing ponazuril can be stored for 28 days with no loss of drug. By applying the FDA Bioanalytical Guidelines (FDA, Guidance for industry: Bioanalytical Method Validation) for method development that states a value should be within 15% deviation from the standard all of the samples were stable for 84 days when stored at -80°C. This study provides much needed information on the stability of ponazuril plasma samples stored at -80°C.

References

- Cox S, Harvill L, Singleton S, Bergman J, DeBolt B. Bioanalytical RP-HPLC method for the determination of ponazuril in plasma. Biomedical Chromatography <https://doi.org/10.1002/bmc.5210>.
- Dirikolu L, Karpiesiuk W, Lehner AF, Tobin T. Toltrazuril sulfone salt: synthesis, analytical detection and pharmacokinetics in the horse. Journal of Veterinary Pharmacology and Therapeutics 2011; 35:265-274.
- FDA, Guidance for industry: Bioanalytical Method Validation 2018. <http://www.fda.gov/downloads/Drugs/GuidanceComplianceRegulatoryInformation/Guidances/UCM070107.pdf>.

Furr M, Kennedy T. Effects of coadministration of corn oil and ponazuril on serum and cerebrospinal fluid concentrations of ponazuril in horses. *Journal of Veterinary Internal Medicine* 2020; 34:1321-1324.

Love D, Gibbons P, Fajt, Jones M. Pharmacokinetics of single-dose oral ponazuril in weanling goats. *Journal of Veterinary Pharmacology and Therapeutics* 2015; 39: 305-308.

Prado ME, Ryman JT, Boileau MJ, Martin-Jimenez T, Meibohm B. Pharmacokinetics of ponazuril after oral administration to healthy llamas (*Lama Glama*). *American Journal of Veterinary Research* 2011; 72:1386-1389.

Stock ML, Elazab ST, Hsu WH. Review of triazine antiprotozoal drugs used in veterinary medicine. *Journal of Veterinary Pharmacology and Therapeutics* 2018; 41:184-194.

Zou M, Guo G, Zhao Y, Zhang Q. Detection, quantifications, and pharmacokinetics of ponazuril in healthy swine. *Journal of Veterinary Pharmacology and Therapeutics* 2014; 37:598-602.

Figure 1

Ponazuril Drug Loss

Percent ponazuril loss in plasma over 84 days (A) all 3 concentrations (B) 0.3 $\mu\text{g/mL}$, (C) 7.5 $\mu\text{g/mL}$, (D) 15 $\mu\text{g/mL}$.




Table 1(on next page)

Ponazuril stability information for 84 days after storage at -80°C.

Concentration results reported in µg/mL; n = 1; * Day 1 samples were not frozen but analyzed immediately after spiking.

1 Table 1. Ponazuril stability information for 84 days after storage at -80°C.

	Concentration at different days of storage												
Starting conc.	Day 1*	Day 7	Day 14	Day 21	Day 28	Day 35	Day 42	Day 49	Day 56	Day 63	Day 70	Day 77	Day 84
0.3	0.3	0.3	0.3	0.3	0.3	0.3	0.291	0.291	0.294	0.294	0.282	0.294	0.273
7.5	7.5	7.5	7.5	7.5	7.5	7.05	7.05	7.2	7.05	7.1	7.0	6.8	6.8
15	15.0	15.0	15.0	15.0	15.0	14.4	14.85	14.85	14.25	14.55	14.7	13.95	14.85
	Percentage change in concentration from starting concentration												
0.3	0	0	0	0	0	0	-3	-3	-2	-2	-6	-2	-9
7.5	0	0	0	0	0	-6	-6	-4	-6	-5	-6	-9	-9
15	0	0	0	0	0	-4	-1	-1	-5	-3	-2	-7	-1

2 Concentration results reported in µg/mL; n = 1; * Day 1 samples were not frozen but analyzed
3 immediately after spiking.

4